

Test tube sealing tape for blood test

118


Features

1. 118 is opacified white tape.
2. It can be stripped and stuck again.
3. It has basic adhesive performance.

Applications

- For sealing of test tubes
- For fixation of sealing materials

Structure



Properties

Product Name	Thickness (mm)	Adhesive force (N/25mm)
118	0.082	7

Precautions for Use

User is responsible for determining whether the KGK product is fit for a particular purpose and suitable for user's method of application. Please remember that many factors can affect the use and performance of a KGK product in a particular application. The materials to be bonded with the product, the surface preparation of those materials, the product selected for use, the conditions in which the product is used, and the time and environmental conditions in which the product is expected to perform are among the many factors that can affect the use and performance of a KGK product. Given the variety of factors that can affect the use and performance of a KGK product, some of which are uniquely within the user's knowledge and control, it is essential that the user evaluate the KGK product to determine whether it is fit for a particular purpose and suitable for the user's method of application. KGK make no warranties on above data.

Precautions for Storage

Please store tapes in a box and cool, dark place, avoiding direct sunlight.
Particularly don't expose to high temperature and humidity.
Over 30°C temperature and 50% or more humidity is strictly prohibited.

KGK Chemical Corporation.
940 Minaminagai Tokorozawa-City saitama-Pref
359-0011 Japan
Tel : +81 4 2944 5151
Mail : info-k@kgk-tape.co.jp
URL : <https://www.kgk-tape.co.jp/>